
Osseous Technologies of America
4500 Campus Drive, Suite 662
Newport Beach, CA 92660 USA
Phone: +1.949.250.0184

Medical Device and QA Services
76 Stockport Road
Timperley, Cheshire, UK WA157SN
Phone: +44 161 711 7181

REF OTA 200 and OTA 205

OSSEOUS SCRAPER INSTRUCTIONS FOR USE

Intended Use: The Osseous Scraper is a manually-operated device which provides a one-step minimally invasive method for harvesting, and collecting cortical and cancellous bone from various sites of the oral cavity and carrying them to the receiving site.

BECAUSE OF RESTRICTIONS WHEN WORKING IN THE POSTERIOR OF THE MOUTH, THE STRAIGHT SCRAPER IS RECOMMENDED FOR THE SYMPHYSIS, RAMUS AND TUBEROSITY AREAS OF THE MOUTH.

THE ANGLED SCRAPER IS DESIGNED TO HARVEST BONE FROM THE EXTERNAL OBLIQUE RIDGE, AND MOLAR AREA.

FIRM PRESSURE IN SHORT (APPROX. 5mm) 'STROKES' OF THE BLADE PRODUCES THE DESIRED RESULTS. **HOLD AT 35 DEGREE ANGLE!**

***NOTE: USE ANY .050 HEX DRIVER AND TWIST THE BLADE EDGE TO EXPOSE NEW CUTTING SURFACE!**

The bone is collected by pushing the end of the device towards the bone surface and simultaneously pulling the device itself backward.

The bone can also be collected with slanting and sideways movements. By varying the angle of the instrument and the pressure exercised, the desired consistency of the bone lamellae or plates can be obtained.

To draw the collected bone, hold the transparent tank facing downward and slide the blade handle backward until completely open.

To transfer the bone lamellae or plates to the receiving site, any type of forceps or spoon can be used.

 **Do not reuse**

 **Sterilized by gamma**

Contents sterile unless package is opened or damaged.

Examine package to ensure it has not been opened and has not been damaged. If the sterile packaging is opened or damaged, do not use and contact Osseous Technologies of America immediately.

After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable laws and regulations.

Explanation of Symbols

	Manufacturer
	European Authorized Representative
	Lot Number
	Catalogue Number
	Do Not Reuse
	Caution, Consult Accompanying Documents
	Sterilization Using Gamma
	Use By
	Quantity

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