


OSSEO-TIP™ BONE GRAFTER

Osseo-Tip™ collects all the bone chips into a removable cup that would otherwise have been wasted during the drilling procedure. After drilling, you simply empty the cup and autogenous bone is ready for immediate use.


The Osseo-Tip™ Bone Grafter proved to be an effective way to harvest sufficient quantities of autogenous bone without the need to employ a bone mill or morselizer. It allowed for bone to be harvested in a safe and rapid manner with no hand fatigue, and without creating a secondary surgical site.

– Dr. Doug-Seok Sohn


Osseo-Tip™ is designed to fit your standard handpiece and is fully Autoclavable.



Place



Drill



Collect